

LAB TOOLS > GLASSWARE

ERLENMEYER

Code: W002



 An Erlenmeyer flask in engineering testing laboratories is used for mixing, heating, and storing liquids or samples. Its conical shape allows for easy swirling without spilling, and it is commonly used in titrations, chemical reactions, and sample preparation, offering both stability and ease of use.

ORDERING INFORMATION

Item	Code
ERLENMEYER - 25 ml	W002G025H
ERLENMEYER - 50 ml	W002G050H
ERLENMEYER - 100 ml	W002G100H
ERLENMEYER - 150 ml	W002G150H





ERLENMEYER - 250 ml	W002G250H
ERLENMEYER - 300 ml	W002G300H
ERLENMEYER - 500 ml	W002G500H
ERLENMEYER - 1000 ml	W002G01KH
ERLENMEYER - 2000 ml	W002G02KH
ERLENMEYER - 3000 ml	W002G03KH
ERLENMEYER - 5000 ml	W002G05KH

OTHER PHOTOS

