



LAB TOOLS > GLASSWARE

DISTILLING FLASK

Code : W004



• A distilling flask in engineering testing laboratories is used for heating liquids to separate components based on their boiling points through distillation. It is commonly used in chemical analysis, purification processes, and the production of solvents or distilled water in the lab.

ORDERING INFORMATION

ltem	Code
DISTILLING FLASK - 100 ml	W004G100H
DISTILLING FLASK - 250 ml	W004G250H
DISTILLING FLASK - 500 ml	W004G500H
DISTILLING FLASK - 1000 ml	W004G01KH