

STIRRING BARS

Code : W005



- Stirring bars in engineering testing laboratories are used for mixing liquids or solutions within containers. Typically made of magnetic material, they are placed inside a flask or beaker and are stirred using a magnetic stirrer. They ensure uniform mixing of samples, chemicals, or reagents during experiments.

ORDERING INFORMATION

Item	Code
STIRRING BARS - Height: 200 mm	W005G200H
STIRRING BARS - Height: 250 mm	W005G250H
STIRRING BARS - Height: 300 mm	W005G300H