

## FILTERING FLASK

Code : W021



- A filtering flask in engineering testing laboratories is used to hold liquids during filtration processes. It is typically used with a vacuum pump to create suction, which speeds up the filtration of liquids through a filter paper, ensuring the separation of solids from liquids in various chemical or material testing procedures.

### **ORDERING INFORMATION**

Item	Code
FILTERING FLASK - 100 ml	W021G100H
FILTERING FLASK - 250 ml	W021G250H
FILTERING FLASK - 500 ml	W021G500H

---

FILTERING FLASK - 1000 ml

W021G01KH

FILTERING FLASK - 2000 ml

W021G02KH

**OTHER PHOTOS**

