

LAB TOOLS > GLASSWARE

PORCELAIN MORTAR and PESTLE

Code: W026



 A porcelain mortar and pestle in engineering testing laboratories are used for grinding, crushing, and mixing solids into fine powders. The non-porous nature of porcelain ensures that no material is absorbed, making it ideal for preparing samples for chemical or material analysis without contamination.

ORDERING INFORMATION

ltem	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H



W026 DATASHEET

PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H

OTHER PHOTOS

