



LAB TOOLS > GLASSWARE

GLASS MORTAR and PESTLE

Code : W028



• A glass mortar and pestle in engineering testing laboratories are used for grinding and mixing materials, particularly for substances that may react with porcelain or other materials. The smooth surface of glass ensures no contamination and is ideal for preparing fine powders or combining chemicals in precise quantities for testing.

ORDERING INFORMATION

| ltem | Code |
|--|-----------|
| PORCELAIN MORTAR and PESTLE - Ø 80 mm | W026C080H |
| PORCELAIN MORTAR and PESTLE - Ø 100 mm | W026C100H |
| PORCELAIN MORTAR and PESTLE - Ø 130 mm | W026C130H |



W028 DATASHEET

| PORCELAIN MORTAR and PESTLE - Ø 160 mm | W026C160H |
|---|-----------|
| PORCELAIN MORTAR and PESTLE - Ø 190 mm | W026C190H |
| PORCELAIN MORTAR and RUBBER-HEADED PESTLE | W027P001H |
| GLASS MORTAR and PESTLE - Ø 100 mm | W028G100H |
| GLASS MORTAR and PESTLE - Ø 120 mm | W028G120H |
| GLASS MORTAR and PESTLE - Ø 150 mm | W028G150H |
| GLASS MORTAR and PESTLE - Ø 180 mm | W028G180H |

OTHER PHOTOS

