


# JAW CRUSHER

Code : M412



- Used to crush samples with medium-hard, hard, brittle and tough materials. Efficiency and safety makes it ideal for sample preparation in laboratories and industrial plants.
- Application examples are so wide such as construction materials, glass, granite, minerals, quartz, rocks, silicon, slag, coal, oxide ceramics ... etc.
- The materials after crushing have high degree of size reduction and high fineness.
- Double action of jaws that opens/closes from the bottom and the crushing mechanism is from top to bottom.
- Supplied complete with electrical motor.

## **TECHNICAL SPECIFICATIONS**

- Double Action of jaws:
  - Opening /closing at bottom

- Crushing from up to bottom
- Jaw opening at the bottom is adjustable.
- 380 V / 50 Hz

### **ORDERING INFORMATION**

| Item                       | Code      |
|----------------------------|-----------|
| JAW CRUSHER (Small)        | M412X001U |
| JAW CRUSHER (Large)        | M412X002V |
| JAW CRUSHER (Extra Large)  | M412X003V |
| SPARE JAWS for M412        | M412P001H |
| SPARE SIDE PLATES for M412 | M412P002H |

### **OTHER PHOTOS**

