

## ■ METHYLENE BLUE SET (AG-041) EN 933-9

- Used to determine the amount of clay in the fine proportions of aggregates.



### ● COMPOSED OF

- 50 ml burette
- Support base with clamp
- Filter paper (Ø 125 mm, pack of 100 pcs)
- Glass Rod (Ø 8 mm x 300 mm)
- 1000 ml beaker
- Methylene blue 100 g
- Kaolinite 500 g
- Mixer with impeller (~ 600 rpm)
  - 220 - 240 V / 50 Hz