

ATTERBERG LIMITS
SOIL
SHRINKAGE LIMIT SET (T-037)

ASTM D427 • AASHTO T-92

- Used to determine the maximum amount of moisture in the sample at which the soil will not shrink while drying.



COMPOSED OF

- | | | |
|-------------------------|-------------------|----------|
| • Shrinkage dish (2 ea) | (Ø 45 mm x 10 mm) | T-037/D |
| • Crystallizing Dish | (Ø 55 mm x 35 mm) | GL-20/S |
| • Shrinkage Prong Plate | | T-037/PP |
| • Evaporating Dish | (Ø 120 mm) | GPP/12 |
| • Graduated Cylinder | (25 ml) | GCM/0025 |
| • Spatula | (100 mm) | GL-17/S |
| • Plastic Carrying Case | | |

LINEAR SHRINKAGE MOULD (T-039)

BS 1377-2

- Used to determine the shrinkage percentage of a soil sample with low clay content.


TECHNICAL SPECIFICATIONS

- Length : 140 mm
- Radius : 12.5 mm