

## **WIRE BASKET**

Code : G100P002H



- Used to hang the aggregate specimens under weighing balance (ET & HT) in Specific Gravity Frame (G-100) to determine the mass of specimen in the air and inside the water to calculate the specific gravity of the specimens
- Consists of stainless-steel frame and wire mesh
- Can be used for various types of tests such as
  - Specific Gravity Determination (Buoyancy Method)
  - Aggregate Soundness Test (Magnesium Sulphate Method)

### **TECHNICAL SPECIFICATIONS**

- Stainless-steel frame and wire mesh

**EQUIPPED WITH**

- Handle

**ORDERING INFORMATION**

Code	Name
G 100 X00X H	SPECIFIC GRAVITY FRAME
G 100 P001 H	CRADLE
G 100 P002 H	DENSITY BASKET (Wire Basket)
G 100 P003 H	WATER TANK
G 100 P004 H	SOUNDNESS BASKET (Mesh Size 1 mm)