

WIRE BASKET

Code : G100P002H



- Used to hang the aggregate specimens under weighing balance (ET & HT) in Specific Gravity Frame (G-100) to determine the mass of specimen in the air and inside the water to calculate the specific gravity of the specimens
- Consists of stainless-steel frame and wire mesh
- Can be used for various types of tests such as
 - Specific Gravity Determination (Buoyancy Method)
 - Aggregate Soundness Test (Magnesium Sulphate Method)

TECHNICAL SPECIFICATIONS

- Stainless-steel frame and wire mesh

EQUIPPED WITH

- Handle

ORDERING INFORMATION

| Code | Name |
|--------------|-----------------------------------|
| G 100 X00X H | SPECIFIC GRAVITY FRAME |
| G 100 P001 H | CRADLE |
| G 100 P002 H | DENSITY BASKET (Wire Basket) |
| G 100 P003 H | WATER TANK |
| G 100 P004 H | SOUNDNESS BASKET (Mesh Size 1 mm) |